

Anti-Sgsh Rabbit Polyclonal Antibody For Wb And Elisa - P51688

Item Number: CM0028373



Introduction

Rabbit polyclonal antibody against human SGSH (N-sulphoglucosamine sulphohydrolase), validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Suitable for lysosomal storage disease research.

[Learn More](#)

Parameter	Value
Product Name	SGSH Rabbit pAb
Remarks/Alias	HSS; SFMD; MP53A; SGSH
Species	Human
Gene ID	6448
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 273-502 of human SGSH (NP_000190.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51688
Protein Weight	57kDa
Shipping	Ice bag