


Anti-Sgta Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O43765

Item Number: CM0029984



Introduction

Anti-SGTA Rabbit Polyclonal Antibody validated for Western blot, IF/ICC, and ELISA applications. Targets the human SGTA co-chaperone involved in protein quality control and tail-anchored protein insertion. Suitable for research use.

[Learn More](#)

Parameter	Value
Product Name	SGTA Rabbit pAb
Remarks/Alias	SGT; Vpu; SGT1; hSGT; alphaSGT; SGTA
Species	Human
GeneID (Human)	6449
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human SGTA (NP_003012.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	O43765
Protein Weight	34kDa
Shipping	Ice bag