


Anti-Sh2D1A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O60880

Item Number: CM0014322



Introduction

SH2D1A Rabbit Polyclonal Antibody validated for Western blot (WB), Immunofluorescence/Immunocytochemistry (IF/ICC), and ELISA. Detects endogenous levels in human and mouse samples. Targets the SLAM-associated protein SAP, a cytoplasmic adapter of 14 kDa. Ideal for NK cell and immune signaling research.

[Learn More](#)

Parameter	Value
Product Name	SH2D1A Rabbit pAb
Remarks/Alias	LYP, SAP, XLP, DSHP, EBVS, IMD5, XLPD, MTCP1, XLPD1, SAP/SH2D1A, SH2D1A
Species	Human
GeneID	4068
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-128 of human SH2D1A (NP_002342.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O60880
Protein Weight	14kDa
Shipping	Ice bag