

Anti-Sh3Bgrl Polyclonal Antibody For Wb, If/Icc, Elisa - O75368

Item Number: CM0022081



Introduction

Rabbit polyclonal antibody against human SH3BGR1 adapter protein. Validated for WB, IF/ICC, ELISA. Cross-reacts with Human, Mouse, Rat. Suitable for studying translation regulation and ubiquitination.

[Learn More](#)

Parameter	Description
Product Name	SH3BGR1 Rabbit pAb
Remarks/Alias	SH3BGR; HEL-S-115; SH3BGR1
Species	Human
GeneID (Human)	6451
GeneID	6451
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-114 of human SH3BGR1 (NP_003013.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75368
Protein Weight	13kDa
Shipping	Ice bag