


Anti-Sh3Bp2 Polyclonal Antibody For Wb, Elisa - P78314

Item Number: CM0017191



Introduction

Rabbit polyclonal antibody targeting human SH3BP2 (3BP2). Validated for Western blot and ELISA, cross-reacts with human and mouse. Recombinant immunogen ensures specific recognition of SH3 domain-binding protein 2 in signal transduction studies.

[Learn More](#)

Parameter	Value
Product Name	SH3BP2 Rabbit pAb
Remarks/Alias	3BP2; CRBM; CRPM; 3BP-2; RES4-23; SH3BP2
Species	Human
GeneID (Human)	6452
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 392-561 of human SH3BP2 (NP_001116153.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P78314
Protein Weight	62kDa
Shipping	Ice bag