

Anti-Sh3GL1 Rabbit Polyclonal Antibody For Wb, Elisa - Q99961

Item Number: CM0027225

Introduction

High-quality rabbit polyclonal antibody targeting human SH3GL1 (Endophilin-A2). Validated for WB and ELISA with cross-reactivity to human, mouse, rat, and monkey. Suitable for endocytosis and membrane trafficking studies.

[Learn More](#)

Parameter	Value
Product Name	SH3GL1 Rabbit pAb
Remarks/Alias	EEN; CNSA1; SH3P8; SH3D2B; SH3GL1
Species	Human
Gene ID (Human)	6455
Gene ID	6455
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 219-368 of human SH3GL1 (NP_003016.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, Monkey
SWISS	Q99961
Protein Weight	41kDa
Shipping	Ice bag