


Anti-Sh3GL2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q99962

Item Number: CM0026142



Introduction

High-quality rabbit polyclonal antibody targeting human SH3GL2 (Endophilin-A1), validated for WB, IF/ICC and ELISA, with mouse cross-reactivity. Ideal for studying clathrin-mediated endocytosis and neuronal development.

[Learn More](#)

Parameter	Value
Product Name	SH3GL2 Rabbit pAb
Remarks/Alias	CNSA2; SH3P4; EEN-B1; SH3D2A; SH3GL2
Species	Human
GeneID (Human)	6456
GeneID	6456
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human SH3GL2 (NP_003017.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q99962
Protein Weight	40kDa
Shipping	Ice bag