


Anti-Sh3GL3 Polyclonal Antibody For Wb, Elisa - Q99963

Item Number: CM0026446



Introduction

Rabbit polyclonal antibody targeting SH3GL3 (Endophilin-A3) for WB and ELISA applications. Cross-reacts with human, mouse, and rat samples. Recombinant full-length protein immunogen ensures specificity.

[Learn More](#)

Parameter	Description
Product Name	SH3GL3 Rabbit pAb
Remarks/Alias	CNSA3; EEN-B2; SH3D2C; SH3P13; HsT19371; SH3GL3
Species	Human
GeneID (Human)	6457
GeneID	6457
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human SH3GL3 (NP_003018.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99963
Protein Weight	39kDa
Shipping	Ice bag