


Anti-Sh3Pxd2A Polyclonal Antibody For Wb, If/Icc, Elisa - Q5Tcz1

Item Number: CM0027011



Introduction

High-quality rabbit polyclonal antibody targeting SH3PXD2A (FISH, TK55). Validated in WB, IF/ICC, ELISA. Cross-reacts with human, mouse. Ideal for invadopodia, podosome, and cancer invasiveness studies.

[Learn More](#)

Parameter	Value
Product Name	SH3PXD2A Rabbit pAb
Remarks/Alias	FISH; TK55; SH3MD1; SH3PXD2A
Species	Human
GeneID (human)	9644
GeneID	9644
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 806-1105 of human SH3PXD2A (NP_055446.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q5TCZ1
Protein Weight	125kDa
Shipping	Ice bag