


Anti-Shank2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Upx8

Item Number: CM0023681



Introduction

High-quality rabbit polyclonal antibody against SHANK2 (SH3 and multiple ankyrin repeat domains protein 2) validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for neuroscience research on synaptic organization and autism spectrum disorders.

[Learn More](#)

Parameter	Value
Product Name	SHANK2 Rabbit pAb
Remarks/Alias	SHANK; AUTS17; CORTBP1; CTTNBP1; ProSAP1; SPANK-3; SHANK2
Species	Human
GeneID (Human)	22941
GeneID	22941
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 500-600 of human SHANK2 (NP_036441.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UPX8
Protein Weight	201kDa
Shipping	Ice bag