

Anti-Shb Rabbit Monoclonal Antibody For Wb, Elisa - Q15464

Item Number: CM0007634

Introduction

Recombinant rabbit monoclonal antibody against human SHB, validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for signal transduction, angiogenesis, and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	SHB Rabbit mAb
Remarks/Alias	bA3J10.2; SHB
Species	Human
GeneID (Human)	6461
GeneID	6461
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 409-509 of human SHB (NP_003019.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15464
Protein Weight	55 kDa
Shipping	Ice bag