

Anti-Shbg Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P04278

Item Number: CM0029574



Introduction

High-quality rabbit polyclonal antibody against human SHBG, validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for hormone transport and receptor-mediated process studies.

[Learn More](#)

Parameter	Value
Product Name	SHBG Rabbit pAb
Remarks/Alias	ABP; SBP; TEBG; SHBG
Species	Human
GeneID (human)	6462
GeneID	6462
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-170 of human SHBG (NP_001031.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04278
Protein Weight	44kDa
Shipping	Ice bag