


Anti-Shc1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P29353

Item Number: CM0030697



Introduction

SHC1 Rabbit pAb is a polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Detects SHC1 protein across human, mouse, and rat samples, targeting the signaling adapter in growth factor and Ras pathways.

[Learn More](#)

Parameter	Value
Product Name	SHC1 Rabbit pAb
Remarks/Alias	SHC; SHCA; SHC1
Species	Human
Gene ID (Human)	6464
Gene ID	6464
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 300-400 of human SHC1 (NP_001189788.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29353
Protein Weight	63kDa
Shipping	Ice bag