

Anti-Shc4 Polyclonal Antibody For Wb, If/Icc, Elisa - Q6S5L8

Item Number: CM0027138

Introduction

High-quality rabbit polyclonal antibody against human SHC4 (RaLP/hShcD). Validated for WB, IF/ICC, and ELISA. Detects human and mouse SHC4, a melanoma-specific adaptor protein involved in tumor migration and invasion.

[Learn More](#)

Specification	Value
Product Name	SHC4 Rabbit pAb
Remarks/Alias	RaLP, SHCD, SHC4
Species	Human
Gene ID (Human)	399694
Gene ID	399694
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-185 of human SHC4 (NP_976224.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q6S5L8
Protein Weight	69kDa
Shipping	Ice bag