

Anti-Shh Rabbit Polyclonal Antibody For Wb, Elisa - Q62226

Item Number: CM0024348



Introduction

Rabbit polyclonal antibody against mouse Sonic hedgehog (SHH) protein, validated for WB and ELISA. Suitable for mouse, human, and rat samples. Ideal for studying developmental signaling and morphogen function.

[Learn More](#)

Parameter	Value
Product Name	SHH Rabbit pAb
Remarks/Alias	Hx; Dsh; Hhg1; Hxl3; ShhNC; M100081; 9530036O11Rik; SHH
Species	Mouse
GeneID	20423
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 187-281 of mouse SHH (NP_033196.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q62226
Protein Weight	48kDa
Shipping	Ice bag