

Anti-Ship1 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q92835

Item Number: CM0007921



Introduction

SHIP1 Rabbit Monoclonal Antibody validated for WB, IF-P, IHC-P, and ELISA. Targets human SHIP1 (Q92835) with cross-reactivity to mouse and rat. Ideal for immune signaling research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	SHIP1 Rabbit mAb
Remarks/Alias	SHIP; SHIP1; SHIP-1; hp51CN; SIP-145; p150Ship
Species	Human
GeneID (Human)	3635
GeneID	3635
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1084-1178 of human SHIP1 (Q92835)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92835
Protein Weight	133kDa
Shipping	Ice bag