

Anti-Shmt2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P34897

Item Number: CM0018067



Introduction

Anti-SHMT2 rabbit polyclonal antibody for WB, IHC-P, ELISA. Detects human, mouse, rat SHMT2. Suitable for studies of one-carbon metabolism, mitochondrial function, and thymidylate biosynthesis.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SHMT2 Rabbit pAb
Remarks/Alias	GLYA; SHMT; mSHMT; NEDCASB; HEL-S-51e; T2
Species	Human
Gene ID (Human)	6472
Gene ID	6472
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 265-504 of human SHMT2 (NP_005403.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P34897
Protein Weight	56kDa
Shipping	Ice bag