


Anti-Shoc2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Uq13

Item Number: CM0028899



Introduction

High-quality rabbit polyclonal antibody against human SHOC2, a leucine-rich repeat scaffold protein in the MAPK pathway. Suitable for western blot, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	SHOC2 Rabbit pAb
Remarks/Alias	SHOC2, SOC2, SUR8, NSLH1, SIAA0862
Species	Human
GeneID (Human)	8036
GeneID	8036
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human SHOC2 (NP_031399.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UQ13
Protein Weight	65kDa
Shipping	Ice bag