

Anti-Shox2 Polyclonal Antibody For Wb, Elisa - O60902

Item Number: CM0026953

Introduction

Rabbit polyclonal antibody against human SHOX2, validated in WB and ELISA, with cross-reactivity to mouse and rat. Suitable for studying developmental regulation and neural somatosensory processing.

[Learn More](#)

Parameter	Value
Product Name	SHOX2 Rabbit pAb
Remarks/Alias	OG12; SHOT; OG12X; SHOX2
Species	Human
GeneID (Human)	6474
GeneID	6474
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60902
Protein Weight	35kDa
Shipping	Ice bag