

Anti-Shp1 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P29350

Item Number: CM0003325



Introduction

High-quality anti-SHP1 rabbit monoclonal antibody validated for WB, IF-P, IHC-P, and ELISA. Detects human, mouse, and rat SHP1 (P29350). Ideal for immunological signaling and hematopoiesis research.

[Learn More](#)

Parameter	Value
Product Name	SHP1 Rabbit mAb
Remarks/Alias	HCP; HCPH; SHP1; SHP-1; HPTP1C; PTP-1C; SHP-1L; SH-PTP1
Species	Human
GeneID (Human)	5777
GeneID	5777
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 496-595 of human SHP1 (P29350).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29350
Protein Weight	68kDa
Shipping	Ice bag