


Anti-Shp1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P29350

Item Number: CM0016687



Introduction

Rabbit anti-SHP1 polyclonal antibody validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat SHP1 (68 kDa; SWISS P29350). Ideal for studying tyrosine phosphatase signaling in immune regulation and hematopoiesis.

[Learn More](#)

Parameter	Value
Product Name	SHP1 Rabbit pAb
Remarks/Alias	HCP; HCPH; SHP1; SHP-1; HPTP1C; PTP-1C; SHP-1L; SH-PTP1
Species	Human
GeneID (Human)	5777
GeneID	5777
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-400 of human SHP1 (NP_002822.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29350
Protein Weight	68kDa
Shipping	Ice bag