


Anti-Shp2 Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q06124

Item Number: CM0016595



Introduction

Anti-SHP2 polyclonal antibody validated for WB, IF/ICC, IP, and ELISA. Cross-reacts with human, mouse, and rat SHP2. Essential for studying MAPK signaling, PD-1/PD-L1 axis, and oncogenic phosphatase function.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SHP2 Rabbit pAb
Remarks/Alias	CFC; NS1; JMML; SHP2; BPTP3; PTP2C; METCDS; PTP-1D; SH-PTP2; SH-PTP3; P2
Species	Human
GeneID (Human)	5781
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-460 of human SHP2 (NP_002825.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q06124
Protein Weight	68kDa
Shipping	Ice bag