


Anti-Shp2 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q06124

Item Number: CM0007369



Introduction

This rabbit monoclonal antibody targets human SHP2 (PTPN11), a key phosphatase in MAPK signaling. Validated for WB, IP, and ELISA, it supports research in signal transduction and disease mechanisms.

[Learn More](#)

Parameter	Value
Product Name	SHP2 Rabbit mAb
Remarks/Alias	CFC; NS1; JMML; SHP2; BPTP3; PTP2C; METCDS; PTP-1D; SH-PTP2; SH-PTP3
Species	Human
Gene ID (Human)	5781
Gene ID	5781
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 502-593 of human SHP2 (Q06124)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	Q06124
Protein Weight	68kDa
Shipping	Ice bag