


Anti-Shp2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q06124

Item Number: CM0025317



Introduction

Rabbit polyclonal antibody targeting human SHP2 (PTPN11), validated for WB, IHC-P, IF/ICC, and ELISA. Recognizes human, mouse, and rat SHP2. Suitable for signal transduction, cancer, and Noonan syndrome research.

[Learn More](#)

Parameter	Value
Product Name	SHP2 Rabbit pAb
Remarks/Alias	CFC; NS1; JMML; SHP2; BPTP3; PTP2C; METCDS; PTP-1D; SH-PTP2; SH-PTP3
Species	Human
GeneID (Human)	5781
GeneID	5781
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 514-593 of human SHP2 (NP_002825.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q06124
Protein Weight	68kDa
Shipping	Ice bag