

Anti-Shp/Nr0B2 Rabbit Polyclonal Antibody For Wb And Elisa - Q15466

Item Number: CM0023404



Introduction

Rabbit polyclonal antibody targeting SHP/NR0B2 (SWISS Q15466), validated for WB and ELISA in human, mouse, and rat. Essential tool for studying nuclear receptor signaling and metabolic circadian rhythm research.

[Learn More](#)

Parameter	Value
Product Name	SHP/NR0B2 Rabbit pAb
Remarks/Alias	SHP; SHP1; SHP/NR0B2
Species	Human
Gene ID	8431
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-257 of human SHP/NR0B2 (NP_068804.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15466
Protein Weight	28kDa
Shipping	Ice bag