

Anti-Siae Rabbit Polyclonal Antibody For Wb And Elisa - Q9Hat2

Item Number: CM0018452

Introduction

High-quality rabbit polyclonal antibody against human SIAE (Sialate O-acetylerase). Suitable for WB and ELISA. Ideal for research on sialic acid metabolism, cell-cell interactions, and tumor biology.

[Learn More](#)

Parameter	Specification
Product Name	SIAE Rabbit pAb
Remarks/Alias	LSE; AIS6; CSEC; YSG2; CSE-C; SIAE
Species	Human
GeneID (Human)	54414
GeneID	54414
Immunogen	Recombinant protein (recombinant fusion protein containing a sequence corresponding to amino acids 270-450 of human SIAE (NP_733746.1)).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9HAT2
Protein Weight	58kDa
Shipping	Ice bag