

Anti-Siah1 Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q8luq4

Item Number: CM0020010



Introduction

Rabbit polyclonal antibody against human SIAH1, validated for WB, IF/ICC, IHC-P, and ELISA. Cross-reacts with mouse and rat. Suitable for ubiquitin-proteasome pathway and cancer research.

[Learn More](#)

Parameter	Value
Product Name	SIAH1 Rabbit pAb
Remarks/Alias	BURHAS; SIAH1A; SIAH1
Species	Human
GeneID (human)	6477
GeneID	6477
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-313 of human SIAH1 (NP_001006611.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IUQ4
Protein Weight	31 kDa
Shipping	Ice bag