

Anti-Siglec-8 Monoclonal Antibody For Fc - Q9Nyz4

Item Number: CM0010045



Introduction

-conjugated rabbit monoclonal antibody against human Siglec-8, validated for flow cytometry. Ideal for detecting Siglec-8 on eosinophils, basophils, and mast cells in allergy and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Siglec-8 mAb
Remarks/Alias	SIGLEC8; SAF2; SIGLEC-8; SIGLEC8L; sialic acid-binding Ig-like lectin 8
Species	Human
Gene ID	27181 (Human)
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 17-363 of human Siglec-8 (NP_055257.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9NYZ4
Protein Weight	44kDa/46kDa/54kDa
Shipping	Ice bag