


Anti-Siglec-8 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nyz4

Item Number: CM0013593



Introduction

Rabbit polyclonal antibody targeting human Siglec-8, an eosinophil and mast cell adhesion receptor. Suitable for WB and ELISA, cross-reactive with mouse and rat. Useful for allergy and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	Siglec-8 Rabbit pAb
Remarks/Alias	SAF2; SIGLEC-8; SIGLEC8L; Siglec-8
Species	Human
Gene ID	27181
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-363 of human Siglec-8 (NP_055257.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NYZ4
Protein Weight	54kDa
Shipping	Ice bag