

Anti-Siglec10 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Lc7

Item Number: CM0022218

Introduction

SIGLEC10 Rabbit Polyclonal Antibody validated for Western blot and ELISA in human and rat samples. Recombinant immunogen, 77kDa protein target. Inquire about bulk quantities and technical consultation from CamelBio.

[Learn More](#)

Parameter	Value
Product Name	SIGLEC10 Rabbit pAb
Remarks/Alias	SLG2; PRO940; SIGLEC-10; SIGLEC10
Species	Human
Gene ID (Human)	89790
Gene ID	89790
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 627-697 of human SIGLEC10 (NP_149121.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q96LC7
Protein Weight	77kDa
Shipping	Ice bag