

Anti-Siglec12 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Pq1

Item Number: CM0029672

Introduction

Rabbit polyclonal antibody against human SIGLEC12 (sialic acid-binding Ig-like lectin 12). Validated for WB and ELISA, cross-reacts with mouse. Recombinant immunogen spanning aa 20-210. For cell adhesion research.

[Learn More](#)

Parameter	Value
Product Name	SIGLEC12 Rabbit pAb
Remarks/Alias	S2V; SLG; SIGLECL1; Siglec-XII; SIGLEC12
Species	Human
GeneID (human)	89858
GeneID	89858
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-210 of human SIGLEC12 (NP_443729.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q96PQ1
Protein Weight	65kDa
Shipping	Ice bag