

Anti-Siglec3/Cd33 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P20138

Item Number: CM0021154



Introduction

Rabbit polyclonal antibody targeting human CD33/SIGLEC3 (Sialic acid-binding Ig-like lectin 3, p67). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. For myeloid cell biology, immune regulation, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	SIGLEC3/CD33 Rabbit pAb
Remarks/Alias	p67; SIGLEC3; SIGLEC-3; SIGLEC3/CD33
Species	Human
GeneID (Human)	945
GeneID	945
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 49-259 of human SIGLEC3/CD33 (NP_001763.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20138
Protein Weight	40kDa
Shipping	Ice bag