


Anti-Siglec5 Rabbit Polyclonal Antibody For Wb And Elisa - O15389

Item Number: CM0030128



Introduction

SIGLEC5 (CD170, OBBP2) rabbit polyclonal antibody validated for Western blot and ELISA, with cross-reactivity to mouse. Used in studying sialic acid-dependent cell adhesion in myeloid lineage cells.

[Learn More](#)

Parameter	Value
Product Name	SIGLEC5 Rabbit pAb
Remarks/Alias	CD170, OBBP2, CD33L2, OB-BP2, SIGLEC-5, SIGLEC5
Species	Human
GeneID (Human)	8778
GeneID	8778
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-434 of human SIGLEC5 (NP_003821.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	O15389
Protein Weight	61kDa
Shipping	Ice bag