

Anti-Siglec6 Rabbit Polyclonal Antibody For Wb, Elisa - O43699

Item Number: CM0018410

Introduction

Validated anti-SIGLEC6 rabbit polyclonal antibody for WB and ELISA. Detects human SIGLEC6 (O43699), a sialic acid-binding lectin mediating cell adhesion, highly expressed in placenta.

[Learn More](#)

Parameter	Specification
Product Name	SIGLEC6 Rabbit pAb
Remarks/Alias	CD327; CD33L; OBBP1; CD33L1; CD33L2; CDW327; SIGLEC6
Species	Human
Gene ID (Human)	946
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-220 of human SIGLEC6 (NP_942142.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O43699
Protein Weight	50kDa
Shipping	Ice bag