

Anti-Siglec9 Polyclonal Antibody For Wb, Elisa - Q9Y336

Item Number: CM0029185

Introduction

Rabbit polyclonal antibody against human SIGLEC9 (CD329). Suitable for Western blot and ELISA. Cross-reacts with human, mouse, and rat. Ideal for studying sialic acid-binding lectin function in immune cells.

[Learn More](#)

Parameter	Value
Product Name	SIGLEC9 Rabbit pAb
Remarks/Alias	CD329; CDw329; FOAP-9; siglec-9; OBBP-LIKE; SIGLEC9
Species	Human
Gene ID (Human)	27180
Gene ID	27180
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 210-340 of human SIGLEC9 (NP_001185487.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y336
Protein Weight	50kDa
Shipping	Ice bag