

Anti-Sik1 Rabbit Polyclonal Antibody For Wb, Elisa - P57059

Item Number: CM0016312

Introduction

SIK1 Rabbit Polyclonal Antibody validated for Western blot and ELISA, detecting human and mouse SIK1. Ideal for studying gluconeogenesis, tumor suppression, circadian rhythms, and muscle differentiation.

[Learn More](#)

Parameter	Value
Product Name	SIK1 Rabbit pAb
Remarks/Alias	MSK; SIK; DEE30; SIK-1; SIK1B; SNF1LK; SIK1
Species	Human
GeneID (Human)	150094
GeneID	150094
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SIK1 (NP_775490.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P57059
Protein Weight	85kDa
Shipping	Ice bag