

## Anti-Sik2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H0K1

Item Number: CM0026109

### Introduction

High-quality rabbit polyclonal antibody against human SIK2 (Q9H0K1), targeting amino acids 827-926. Validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for metabolic and autophagy research.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | SIK2 Rabbit pAb                                                                                                   |
| Remarks/Alias    | QIK; SIK-2; SNF1LK2; LOH11CR1I; SIK2                                                                              |
| Species          | Human                                                                                                             |
| Gene ID (Human)  | 23235                                                                                                             |
| Gene ID          | 23235                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 827-926 of human SIK2 (NP_056006.1) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, ELISA                                                                                                         |
| Cross Reactivity | Human, Mouse, Rat                                                                                                 |
| SWISS            | Q9H0K1                                                                                                            |
| Protein Weight   | 104kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                           |