

Anti-Sik3 Polyclonal Antibody For Wb And Elisa - Q9Y2K2

Item Number: CM0028952



Introduction

Rabbit anti-SIK3 polyclonal antibody validated for Western blot and ELISA in human, mouse, and rat. Targets serine/threonine-protein kinase SIK3 involved in mTOR signaling and chondrocyte differentiation. Suitable for kinase research.

[Learn More](#)

Parameter	Value
Product Name	SIK3 Rabbit pAb
Remarks/Alias	L19, QSK, SEMDK, SIK-3, SIK3
Species	Human
GeneID (Human)	23387
GeneID	23387
Immunogen	Recombinant protein corresponding to amino acids 60-310 of human SIK3 (NP_079440.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y2K2
Protein Weight	145kDa
Shipping	Ice bag