

Anti-Sin3A Monoclonal Antibody For Wb, Elisa - Q96St3

Item Number: CM0006919

Introduction

Rabbit monoclonal antibody against human Sin3A (SWISS: Q96ST3) for WB and ELISA applications. Validated cross-reactivity with human and mouse. Ideal for transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	Sin3A Rabbit mAb
Remarks/Alias	WITKOS; DEL15Q24; CHR15DELq24; Sin3A
Species	Human
Gene ID	25942
Immunogen	Synthetic peptide (a synthetic peptide corresponding to a sequence within amino acids 1-100 of human Sin3A (NP_056292.1))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96ST3
Protein Weight	145kDa
Shipping	Ice bag