

Anti-Sin3A Polyclonal Antibody For Wb, Elisa - Q96St3

Item Number: CM0013276

Introduction

Rabbit anti-Sin3A polyclonal antibody validated for WB and ELISA; cross-reacts with human and rat Sin3A (~145 kDa). Suitable for chromatin, circadian rhythm, and neurodevelopmental studies.

[Learn More](#)

Parameter	Value
Product Name	Sin3A Rabbit pAb
Remarks/Alias	WITKOS; DEL15Q24; CHR15DELq24; Sin3A
Species	Human
GeneID (Human)	25942
GeneID	25942
Immunogen	Recombinant protein corresponding to amino acids 1-120 of human Sin3A (NP_056292.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q96ST3
Protein Weight	145kDa
Shipping	Ice bag