

## Anti-Sin3A Rabbit Polyclonal Antibody For Wb - Q96St3

Item Number: CM0019589

### Introduction

Rabbit polyclonal antibody against human Sin3A transcriptional corepressor, validated for WB, IF/ICC, IP, ELISA. Detects human, mouse, rat. Suitable for studying transcriptional repression and circadian rhythms.

[Learn More](#)

| Parameter        | Value                                                                                                                    |
|------------------|--------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Sin3A Rabbit pAb                                                                                                         |
| Remarks/Alias    | WITKOS; DEL15Q24; CHR15DELq24; Sin3A                                                                                     |
| Species          | Human                                                                                                                    |
| Gene ID (Human)  | 25942                                                                                                                    |
| Gene ID          | 25942                                                                                                                    |
| Immunogen        | Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-120 of human Sin3A (NP_056292.1). |
| Source           | Rabbit                                                                                                                   |
| Category         | Polyclonal Antibodies                                                                                                    |
| Application      | WB, IF/ICC, IP, ELISA                                                                                                    |
| Cross Reactivity | Human, Mouse, Rat                                                                                                        |
| SWISS            | Q96ST3                                                                                                                   |
| Protein Weight   | 145kDa                                                                                                                   |
| Shipping         | Ice bag                                                                                                                  |