

## Anti-Sipa1 Rabbit Polyclonal Antibody For Wb And Elisa - Q96Fs4

Item Number: CM0016221

### Introduction

Anti-SIPA1 rabbit polyclonal antibody for western blot and ELISA applications. Detects human and mouse SIPA1 (p130 SPA-1, 112 kDa, Uniprot Q96FS4), immunogen amino acids 763-1042.

[Learn More](#)

| Parameter        | Value                                                                                                                 |
|------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name     | SIPA1 Rabbit pAb                                                                                                      |
| Remarks/Alias    | SPA1; SIPA1                                                                                                           |
| Species          | Human                                                                                                                 |
| Gene ID (Human)  | 6494                                                                                                                  |
| Gene ID          | 6494                                                                                                                  |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 763-1042 of human SIPA1 (NP_694985.29). |
| Source           | Rabbit                                                                                                                |
| Category         | Polyclonal Antibodies                                                                                                 |
| Application      | WB, ELISA                                                                                                             |
| Cross Reactivity | Human, Mouse                                                                                                          |
| SWISS            | Q96FS4                                                                                                                |
| Protein Weight   | 112 kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                               |