

Anti-Sirt1 Polyclonal Antibody For Wb, Elisa - Q923E4

Item Number: CM0028493

Introduction

Rabbit polyclonal antibody against mouse SIRT1 (Sir2), validated for WB and ELISA. NAD-dependent deacetylase involved in metabolism, aging, apoptosis, and circadian rhythms. Ideal for epigenetic and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	SIRT1 Rabbit pAb
Remarks/Alias	Sir2; Sir2a; SIR2L1; Sir2alpha
Species	Mouse
GeneID (human)	23411
GeneID	93759
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 638-737 of mouse SIRT1 (NP_062786.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q923E4
Protein Weight	80 kDa/60 kDa
Shipping	Ice bag