

Anti-Sirt1 Polyclonal Antibody For Wb, Ip, Ihc-P, Elisa - Q96Eb6

Item Number: CM0014944



Introduction

High-quality rabbit polyclonal antibody against human SIRT1, validated for WB, IP, IHC-P, and ELISA. Suitable for studying apoptosis, metabolism, and chromatin regulation. Cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	SIRT1 Rabbit pAb
Remarks/Alias	SIR2; SIR2L1; SIR2alpha; SIRT1
Species	Human
GeneID (Human)	23411
GeneID	23411
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 448-747 of human SIRT1 (NP_036370.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96EB6
Protein Weight	82 kDa
Shipping	Ice bag