

Anti-Sirt1 Rabbit Monoclonal Antibody (Ko Validated) For Wb, If/Icc, Elisa - Q96Eb6

Item Number: CM0003213



Introduction

High-quality rabbit monoclonal antibody against human SIRT1, validated for WB, IF/ICC, and ELISA. Ideal for studying NAD-dependent deacetylation, apoptosis, metabolism, DNA repair, and gene regulation.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SIRT1 Rabbit mAb
Remarks/Alias	SIR2; SIR2L1; SIR2alpha; T1
Species	Human
Gene ID (Human)	6772
Gene ID	23411
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q96EB6
Protein Weight	82kDa
Shipping	Ice bag