

Anti-Sirt2 Rabbit Monoclonal Antibody For Wb, Elisa - Q8Ixx6

Item Number: CM0009545



Introduction

High-quality rabbit recombinant monoclonal antibody targeting human SIRT2 (NAD-dependent deacetylase sirtuin-2). Validated for Western blot (WB) and ELISA. Ideal for studies on cell cycle, epigenetics, and tumor suppression.

[Learn More](#)

Parameter	Value
Product Name	SIRT2 Rabbit mAb
Remarks/Alias	SIR2; SIR2L; SIR2L2; SIRT2
Species	Human
GeneID	22933
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SIRT2 (Q8IXJ6).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8IXJ6
Protein Weight	43kDa
Shipping	Ice bag