

Anti-Sirt2 Rabbit Polyclonal Antibody (Pab) For Wb, Ihc-P, If/Icc, Elisa - Q8lxj6

Item Number: CM0017735



Introduction

SIRT2 Rabbit pAb validated for WB, IHC-P, IF/ICC, ELISA. Targets human SIRT2 (NAD-dependent deacetylase) involved in cell cycle, autophagy, and microtubule dynamics. Cross-reacts with mouse, rat.

[Learn More](#)

Parameter	Value
Product Name	SIRT2 Rabbit pAb
Remarks/Alias	SIR2; SIR2L; SIR2L2; SIRT2
Species	Human
Gene ID (Human)	22933
Gene ID	22933
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-350 of human SIRT2 (NP_036369.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8LXJ6
Protein Weight	43kDa
Shipping	Ice bag