

Anti-Sirt3 Rabbit Monoclonal Antibody (Ko Validated) For Wb, Ihc-P, Ip, Elisa - Q9Ntg7

Item Number: CM0003958



Introduction

High-specificity anti-SIRT3 rabbit monoclonal antibody, KO validated, for WB, IHC-P, IP, ELISA. Targets human and rat SIRT3, a mitochondrial NAD-dependent protein deacetylase regulating energy metabolism.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SIRT3 Rabbit mAb
Remarks/Alias	SIR2L3; T3
Species	Human
Gene ID (Human)	23410
Gene ID	23410
Immunogen	Synthetic peptide; A synthetic peptide corresponding to a sequence within amino acids 300-399 of human SIRT3 (NP_036371.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9NTG7
Protein Weight	44kDa
Shipping	Ice bag