

Anti-Sirt3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q9Ntg7

Item Number: CM0029761



Introduction

Rabbit polyclonal antibody against human SIRT3 (SIR2L3), validated for WB, IHC-P, IF/ICC, IP, ELISA. Detects human, mouse, rat SIRT3, a mitochondrial NAD-dependent deacetylase involved in metabolism and stress response. Ideal for energy metabolism research.

[Learn More](#)

Parameter	Value
Product Name	SIRT3 Rabbit pAb
Remarks/Alias	SIR2L3; SIRT3
Species	Human
Gene ID	23410 (Human)
Immunogen	Recombinant protein A synthetic peptide corresponding to a sequence within amino acids 388-399 of human SIRT3 (NP_036371.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NTG7
Protein Weight	44kDa
Shipping	Ice bag