

Anti-Sirt3 Rabbit Polyclonal Antibody [Ko Validated] For Wb And Elisa - Q9Ntg7

Item Number: CM0029154



Introduction

[KO Validated] SIRT3 rabbit polyclonal antibody for Western blot and ELISA applications. Validated in human, mouse, and rat. Suitable for studying mitochondrial deacetylase function and energy metabolism.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SIRT3 Rabbit pAb
Remarks/Alias	SIR2L3; T3
Species	Human
GeneID	23410
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 299-399 of human SIRT3 (NP_036371.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NTG7
Protein Weight	44kDa
Shipping	Ice bag